

# NurseAdvocate™ AIO IS 100% FIELD SERVICEABLE

- Drastically reduce downtime of the mobile workstation away from the floors
- Very little interruption of Clinician workflow
- Minimal IT staff labor hours are required to service the NurseAdvocate™ AIO computer
- Your facility saves time and money, avoiding unnecessary shipping fees by servicing the AIO on-site

---

## NurseAdvocate™ AIO IS A SMART INVESTMENT



### Easy Access To All Internal Components:

- Two security screws are the only hardware to remove for quick, easy access to service the NurseAdvocate™ AIO
- Takes less than 1 minute to remove both screws with a simple Allen wrench



### Unique Hinged Design:

- Computer casing opens within seconds for easy access to all internal components
- NurseAdvocate™ AIO does NOT need to be removed from the mobile workstation
- NurseAdvocate™ AIO can be serviced anywhere in your facility



### Industry Standard, Readily Available Componentry:

- Non-proprietary components significantly reduce cost
- Stinger Medical will provide your facility with a list of specific parts and part numbers for easy ordering as necessary
- We recommend your IT support staff maintains a small inventory of components to have ready for immediate use
- Rapid turnaround time on device support ensures the mobile workstation is back on the floor ASAP for Clinicians

---

### SERVICE SOLUTIONS

- Maintenance Programs
- Integration Services
- Imaging Services
- Extended Warranties

### MADE ESPECIALLY FOR

- Healthcare Environments
- Mobile Clinical Workstations
- Fixed Wall Mounts
- Fixed Ceiling Mounts

### POWER SOLUTIONS

- MOBIUS POWER™
- NurseAlert™ Power System